


# ergo

## File Caddy ...

Autoclavable foam for storage and cleaning of rotary files while in motion.

Autoclavable foam for storage and cleaning of hand files.



Load up to 5 rotary files for a quick change directly onto endo motor.

Two rulers to quickly reference working length.

Autoclavable bracelet is comfortable to wear and can be adjusted to size.

## FOR THE ORGANIZATION AND CLEANING OF ENDODONTIC FILES.

Designed in collaboration with the experienced team at Style Italiano Endodontics.

Ergo helps improve efficiency and ergonomics of routine endodontic therapy.

### **Ergo - [800000] - \$75.00**

- (1) Bracelet
- (3) Gray Foam Inserts
- (3) Blue Foam Inserts

### **Foam Insert 48/pk - [800001] - \$40.00**

- (24) Gray Foam Inserts
- (24) Blue Foam Inserts